

Deltalab Sterilisation Indicator Tape (Ethylene Oxide)

Product Code: I91049

Brand: Deltalab

Deltalab's Sterilisation Indicator Tape (Ethylene Oxide) is self-adhesive and made of semi-creped paper. It also sticks on different surfaces such as plastic, metal, fabrics, paper, cardboard, etc. Will remain affixed for 18 months.

Features:

- Overall thickness of 0.11 +/- 0.1 mm
- Indicator for sterilisation by ethylene oxide gas. The word GAS printed in brown changes to green once it has been sterilised as per the following specifications: exposition to 600 mg/l EO and 60%RH for minimum 20 minutes at 54 °C, and exposition to 600 mg/l EO and 60%RH for minimum 25 minutes at 37 °C

Specifications:

Type of Sterilisation	Ethylene Oxide Gas
Length (mm)	50
Width (mm)	19
Sterility	Non-Sterile
Quantity	1

