

Medline 1.5ml Microcentrifuge Tubes

Product Code: 200400SP

Brand: Medline Scientific

Made from selected high-quality and biologically inert polypropylene, Medline's 1.5ml Microcentrifuge Tubes are suitable for use at 121°C/2.1 bar (cap not closed).

The leak-free caps on these graduated microcentrifuge tubes can be sealed efficiently and will stay closed during centrifugation with excellent visibility.

Come in packs of 1000.

Features:

- Made from selected high-quality and biologically inert polypropylene.
- High resistance to temperatures and chemicals.
- Autoclavable at 121°C/2.1 bar (cap not closed).
- Place for markings on label area and flat cover.
- Perfect closure of leak free cap due to the design of the lid.
- Membrane can be punctured by sampling needles.
- Graduation.
- Round inside bottom for easy recovery of pellets or fluids.
- Easy closure (e.g. for use in chemical laboratories).

