

## Deltalab Sterile Amies Charcoal Swabs (Snappable Polystyrene & Viscose)

**Product Code:** 300285

**Brand:** Deltalab

Deltalab's Sterile Amies Charcoal Swabs (Snappable Polystyrene & Viscose) are a modification of Cary Blair media, which is at the same time a modification of Stuart media. It assures the viability of organisms such as *Trichomonas Vaginalis*, *Streptococcus Pyogenes*, and *Streptococcus Pneumoniae*.

Swabs are sterilised by radiation.



### Specifications

<b>Peel-Pack Dimensions (mm)</b>	38 x 210
<b>Material(s)</b>	Snappable Polystyrene and Viscose (with charcoal)
<b>Expiry Date</b>	30 months after sterilisation date
<b>Quantity</b>	100