

Deltalab Graduated Beakers (100ml)

Product Code: 409100

Brand: Deltalab

Ideal for general laboratory use, Deltalab's Graduated Beakers (100ml) are made of durable polypropylene for safe use. These lab beakers are autoclavable up to 121°C, have clearly printed graduations, and are in compliance with ISO 6706 and BS 5404 standards.

Specifications

Capacity (ml)	100
Diameter (mm)	53
Height (mm)	72
Material	Polypropylene
Sterility	Non-Sterile
Quantity	12

