



LSFSF107UK 107L Laboratory Essential Freezer – Solid

LSFSF107UK

- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms

SPECIFICATION

| | |
|----------------------------------|-----------------|
| Dimensions | 495 x 620 x 830 |
| Weight | 33kg |
| Gross Capacity (Litres) | 107 Litres |
| Net Capacity (Litres) | 103 Litres |
| Unit External Width (mm) | 495mm |
| Unit External Depth (mm) | 620mm |
| Unit External Height (mm) | 830mm |
| Product Weight (kg) | 33kg |
| Packed Width (mm) | 560mm |
| Packed Depth (mm) | 850mm |
| Packed Height (mm) | 655mm |
| Packed Weight (kg) | 36.5kg |
| Country Of Origin | UK |
| Internal Dimensions (mm) (WxDxH) | 430 x 700 x 440 |
| Temperature Range | -18°C to -25°C |
| Colour | White/Green |

| | |
|------------------------------|-------------------------|
| External Finish | Painted Steel |
| Internal Finish | HIPS* |
| Type Of Plastic | SRL800 (HIPS/ABS Mix) |
| Door Finish | Painted Steel |
| Door Type | Solid |
| Door Style | Hinged Reversible |
| Door Liner | Flat |
| Door Open Depth (mm) | 1050mm |
| Replacement Door Seal | No |
| Open Door Sensor | Yes |
| Lock | Yes (2 Keys) |
| Adjustable Feet | 1 Adjustable Front Foot |
| Porthole | No |
| Temperature Display On Panel | Yes |
| Temperature Alarm | Yes |
| Remote Alarm Contacts | No |
| Audible Alarms | Yes |
| Visual Alarms | Yes |
| Number Of Shelves | 2 |
| Shelf Size (mm) (WxDxH) | 427 x 374 x 6 |
| Shelf Material | Polythene Coated Wire |
| Number Of Baskets | 1 |
| Basket Material | Polythene Coated Wire |
| Basket Size (mm) (WxDxH) | 395 x 180 x 150 |
| Interior Light | No |
| Fan | No |
| Spark Free | Yes |
| Defrost Setting | Manual |
| Controller Type | Digital External |
| Number Of Probes | 1 |

| | |
|----------------------------|--------------------------------|
| Probe Type | Air |
| Ambient Temperature Range | 10°C to 32°C |
| Refrigerant Type | R600a |
| CFC Free | Yes |
| HFC Free | Yes |
| HCFC Free | Yes |
| Ammonia Free | Yes |
| UK 3 Pin Plug | Yes |
| Other Plug Types Available | Yes |
| Voltage Supply (V) | 230V |
| Insulation Thickness (mm) | 35mm |
| Insulation Type | High Density Polyurethane Foam |
| Blowing Agent | Cyclopentane |
| Anti Microbial Handle | Yes |
| Wall Brackets | No |
| Optional Drawer System | No |
| Optional Digital Lock | No |
| Other Door Option | No |



For more information contact us on
 email: enquiries@medlinescientific.com