

Medline Analogue Hotplate (MS-100)

Product Code: MS100

Brand: Medline Scientific

Medline's Analogue Hotplate (MS-100) offers a large heating area of 180 x 180mm with safe controlled heating up to 400°C. Constructed from chemically resistant ceramic coated steel, this lab hotplate delivers excellent reliability with reproducible results at a sensible price. A standard 12.5mm screwed thread is included for the addition of a retort rod (not included).



Specifications

Description	Hotplate
Plate Area (mm)	180 x 180
Temperature Range (°C)	Max. 400
Dimensions (W x D x H)	410 x 210 x 150mm
Weight (kg)	4