

## Medline Polystyrene Blue Dipper Samplers with Handles (100ml)

**Product Code:** 309688

**Brand:** Medline Scientific

Medline's Polystyrene Blue Dipper Samplers with Handles (100ml) allow laboratories to carry out simple sample collections and subsequent transportation. The handle neatly snaps off after sample collection from the vat for the ease of transportation.

Each dipper is supplied with a flow-seal screw cap to ensure leak-proof closure.

### Specifications

<b>Colour</b>	Blue
<b>Volume (ml)</b>	100
<b>Handle Length (mm)</b>	330
<b>Material</b>	Polystyrene
<b>Inner Pack</b>	1
<b>Quantity</b>	50

