

Medline Metal Cap Containers with Printed Labels (150ml)

Product Code: 309701M

Brand: Medline Scientific

Medline's Metal Cap Containers with Printed Labels (150ml) are sterile containers with a metal flowed seal inert liner cap.

Tray wrapped, these effective sample containers are manufactured from medical grade virgin polystyrene and the caps guarantee samples to be safely secured with their leak proof seal.

Specifications:

Volume (ml)	150
Labels	Printed
Material	Polystyrene
Sterility	Sterile
Quantity	120

